

# PURECLEANSE

# **Instructions For Use**

# **Dispensing Cap**

BIHOCL PureCleanse is a safe irrigation, moisturizing, and debridement solution formulated to actively cleanse wounds and prepare the wound bed for healing in a safe, effective and natural way.



The BIHOCL PureCleanse dispensing cap is a versatile applicator with many use cases.



#### Cleanse the wound

Apply BIHOCL PureCleanse liberally to irrigate wound and peri-wound tissue. Ensure the wound is thoroughly saturated. No rinse is required.



Soak the wound

Saturate gauze with BIHOCL PureCleanse solution and apply to wound, soaking for ~5 minutes. Remove gauze and use it to wipe wound bed for further debridement.



Pack the wound

Saturate packing gauze with BIHOCL PureCleanse solution and place it into the wound as appropriate.

### Indications

BIHOCL PureCleanse is indicated for the management of wounds such as Stage I-IV pressure ulcers, partial and full-thickness wounds, diabetic foot ulcers, post-surgical wounds, 1st/2nd degree burns, grafted and donor sites, and minor irritations of the skin in addition to moistening and lubricating absorbent wound dressings. The product is safe for use on adult and adolescent patients.

### Contraindications

Avoid use on individuals allergic or sensitive to hypochlorous acid (HOCI).

### Storage

Store BIHOCL PureCleanse in a room temperature environment and away from direct sunlight, heat, or extreme cold. Do not freeze. For adequate efficacy, do not use the solution past its expiration date (printed on the bottle itself).

## Disposal

Dispose of components such as the BIHOCL probe, shield, and/or LiquidLock super absorbent pad into a biohazard waste container. Dispose of remaining BIHOCL PureCleanse solution and plastic bottle in the regular garbage disposal, and please recycle wherever possible.



# BIHOCL

# **PURECLEANSE**

# **Instructions For Use**



# Probe / Shield



### Prepare the wound

Assess the wound for signs of infection, and evaluate the wound bed for granulation tissue, slough, and exudate. Remove as much as possible prior to treatment with BIHOCL PureCleanse.



#### Cleanse the wound

Manually compress the bottle to produce a direct, consistent stream of solution onto the wound if shield is used, or into the wound if probe is used.



### Prepare the product

Open the packaging of the desired attachment for treatment (probe, shield, or both). Firmly screw the attachment onto the bottle of BIHOCL PureCleanse solution.



### Dispose responsibly

Once treatment is complete, refer to disposal instructions to safely discard all delivery system components and any excess solution.



#### Position the bottle

When ready to begin cleansing, position the bottle over the wound and invert the bottle. Place the shield or probe as close as possible to the wound, ensuring complete coverage and minimal splash-back.

When to use **PURECLEANSE** 



Questions, comments, or concerns?