



Pure Hypochlorous Acid Cleanser for Wounds, Skin & Burns

BIHOCL PureCleanse is intended for cleansing, irrigating, moistening, debridement, and removal of foreign material including microorganisms and debris from exuding and/or dirty wounds, acute and chronic dermal lesions such as stage I-V pressure ulcers, stasis ulcers, diabetic ulcers, post-surgical wounds, first and second degree burns, abrasions, minor irritations of the skin, diabetic foot ulcers, ingrown toenails, grafted and donor sites, and exit sites. It is also intended for moistening and lubricating absorbent wound dressings.

Available Now
on All National GPO Contracts



Contract Award #06193
AWC - Wound Cleanser & Debrider, Sodium Hypochlorite Acid or Hypochlorous Acid



Contract Award #CW16063 (Mohawk Medbuy)
Contract Award #C11825 (Plexxus)
Wound Care



DISPENSING CAP

A versatile applicator for cleansing, soaking, or packing a wound with saturated gauze.

REF	PRODUCT DESCRIPTION	VOLUME	EA / BOX	HEALTHPRO CONTRACT #	MOHAWK MEDBUY / PLEXUS CONTRACT #
33055-118DC	Hypochlorous Acid Wound Cleanser 330 PPM 5.5 pH Bottle w/ Dispensing Cap Non-Latex Non-Sterile Single Patient Use Disposable	118 mL	24	06193	CW16063 / C11825
33055-250DC	Hypochlorous Acid Wound Cleanser 330 PPM 5.5 pH Bottle w/ Dispensing Cap Non-Latex Non-Sterile Single Patient Use Disposable	250 mL	24	06193	CW16063 / C11825
33055-475DC	Hypochlorous Acid Wound Cleanser 330 PPM 5.5 pH Bottle w/ Dispensing Cap Non-Latex Non-Sterile Single Patient Use Disposable	475 mL	16	06193	CW16063 / C11825
33055-1000DC	Hypochlorous Acid Wound Cleanser 330 PPM 5.5 pH Bottle w/ Dispensing Cap Non-Latex Non-Sterile Single Patient Use Disposable	1000 mL	9	06193	CW16063 / C11825



TRIGGER SPRAYER

An efficient, easy-to-use applicator for in-clinic treatment or at-home patient use.

REF	PRODUCT DESCRIPTION	VOLUME	EA / BOX	HEALTHPRO CONTRACT #	MOHAWK MEDBUY / PLEXUS CONTRACT #
33055-118TS	Hypochlorous Acid Wound Cleanser 330 PPM 5.5 pH Bottle w/ Trigger Sprayer Non-Latex Non-Sterile Single Patient Use Disposable	118 mL	24	06193	CW16063 / C11825



SPIKEABLE INSTILLATION CAP

A high quality silicone seal designed to be easily spiked while preventing leakage during use.

REF	PRODUCT DESCRIPTION	VOLUME	EA / BOX	HEALTHPRO CONTRACT #	MOHAWK MEDBUY / PLEXUS CONTRACT #
33055-250SC	Hypochlorous Acid Wound Cleanser 330 PPM 5.5 pH Bottle w/ Spikeable Cap Non-Latex Non-Sterile Single Patient Use Disposable	250 mL	24	06193	CW16063 / C11825
33055-475SC	Hypochlorous Acid Wound Cleanser 330 PPM 5.5 pH Bottle w/ Spikeable Cap Non-Latex Non-Sterile Single Patient Use Disposable	475 mL	16	06193	CW16063 / C11825
33055-1000SC	Hypochlorous Acid Wound Cleanser 330 PPM 5.5 pH Bottle w/ Spikeable Cap Non-Latex Non-Sterile Single Patient Use Disposable	1000 mL	9	06193	CW16063 / C11825



We're here to help!

sales@biomiq.health | 1-800-573-1295 | www.Biomiq.health/WoundCare

BIOMIQ