

PURECLEANSE™ STAT

FOR TRAUMA & URGENT CARE

REF STAT250

HEALTH CANADA MEDICAL DEVICE
LICENSE NO. 111117

Rapid Antimicrobial Cleansing Lavage & Debridement System for Wounds, Skin and Burns

Shield your patients from the hidden threat of harmful microorganisms. In the often urgent, high-pressure environments where you provide care, these invisible invaders lurk everywhere—on surfaces, skin, and within the body's natural flora. Though unseen, their impact can be severe, making the right choice of wound irrigation solution crucial for optimal patient care.

Flushing the area with pure hypochlorous acid (HOCl) effectively removes debris, contaminants, and bacteria, reducing the risk of infection and kickstarting the healing process.

PureCleanse™ STAT is the ideal choice for fast, effective antimicrobial cleansing and irrigation debridement for traumatic or urgent wounds, skin complications and acute burns.

It's as effective as it is safe.

Proudly manufactured in Canada.



When your patient needs fast, effective and safe wound irrigation, you need PureCleanse™ STAT.

QuickTab™ Technology



**Twist.
Cleanse.
Done.**



Eazy Squeeze

Ergonomic Grip



Why choose PureCleanse™ STAT?

Broad-spectrum antimicrobial wound cleansing
Non-toxic, non-sensitizing, non-irritating, and fully tissue compatible
Optimal pressure jet lavage system (up to 12 PSI)
A no-rinse, no-sting solution with the same safety profile as saline
Supported by independent in-vivo and in-vitro studies

Ideal for use in:

Urgent Care Clinics
Emergency & Trauma Rooms
Wound Clinics
First Responder Care
Military Medical Units

Indicated for:

Partial and full-thickness wounds
Skin abrasions
1st/2nd Degree Burns
Post-surgical wounds
Acute dermal lesions



Add PureCleanse™ STAT to your
urgent care irrigation plan today!

sales@biomiq.health
1 800 573 1295
www.biomiq.health

BIOMIQ

Manufactured for Biomiq Inc. by Sterasure Inc.
137 Glasgow St., Unit 115, Kitchener,
Ontario N2G 4X8, Canada
Rev001

PURECLEANSE™

STAT

Cleanse In Any Crisis.

Backed by Third-Party Laboratory Testing

PureCleanse STAT™ is unsurpassed in its combination of safety, efficacy, and stability.



Formulation and quality tested and validated in numerous studies by the University of Waterloo.*



Purity tested by the National Institutes of Health (NIH).*

Quick and Effective for Optimal Infection Protection

In-Vitro Time Kill Assay Test Results Measuring Pathogenic Colony Log**
Reductions of BIHOCL PureCleanse™ Hypochlorous Acid (HOCl) Solution

Pathogen	Elimination
Acinetobacter baumannii	>99.999%
Aspergillus niger spores	99.9999%
Biofilm, polymicrobial	99.96%
Candida albicans	99.9999%
Candida auris	99.99997%
Coronavirus SARS-CoV-2	99.994%
Enterobacter cloacae	>99.9999%
Escherichia coli NDM-1	>99.99999%
Human papillomavirus (HPV 16 & 18)	>99.999%
Poliovirus	99.9999%
Prions (vCJD, others)	>99.9999%
Pseudomonas aeruginosa	99.9999%
Staphylococcus aureus, methicillin-resistant (MRSA)	>99.999%

*Data on file.

**This is a non-exhaustive list of pathogens tested.
Full data available upon request.

How does PureCleanse™ STAT compare? The choice is clear.

	PURECLEANSE STAT	Normal Saline	Sodium Hypochlorite (Dakin's Solution)	Chlorhexidine Gluconate (CHG)	Hydrogen Peroxide	Povidone Iodine	PHMB	Silver Nitrate
Kills Full Range of Bacteria, Viruses, Spores, and Fungi	●		●			●		
Non-Cytotoxic	●	●						
No Rinse Required	●	●						●
Safe for Mucosal Tissue	●	●				●	●	
Disrupts Biofilm	●		●	●	●		●	●
Biomimetic Solution	●							
Environmentally Friendly	●	●						

How to Order:

Available Now on National GPO Contracts



Order Through Your Preferred Distribution Partner

Contact your preferred distributor
to place an order, or visit our
website for more information.

Format

250ml / EA
20 EA / BX

Distributor

Biomiq (Direct)

Cardinal Health

Medline Canada

The Stevens Company

Ordering Code

STAT250

STAT250

G28-STAT250

924-STAT250